


Alexer[®]

Protective backing layer

Breathable and comfortable non-woven helps prevent leakage



Distribution layers

Hold and transmit the exudates into absorbent layer

Super absorbent core

Absorbs and holds up to 20times its own weight

Wound contact layer

Soft and comfortable non-adherent surface minimized trauma and pain during dressing changes

SUPER ABSORBENT DRESSING


◇ Non Adherent



DESCRIPTION

Alexer® Super absorbent dressing-non adherent is indicated for management of moderately to heavily exuding wounds, including: leg ulcers, pressure ulcers, diabetic foot ulcers and dehisced surgical wounds.


APPLICATIONS




Super absorbent polymer core absorbs and locks up to 20 times its own weight.

Minimize the risk of maceration and leakage.

Non-adherent Conform to body contours.



Type	Size	Package	Ref.No.
 Square	10X10cm	5pcs/bx,40bxes/ctn	S50021001
	7.5X7.5cm	5pcs/bx,40bxes/ctn	S50011002
	12.5X12.5cm	5pcs/bx,40bxes/ctn	S50021003
	15X15cm	5pcs/bx,40bxes/ctn	S50011004
	10X20cm	5pcs/bx,40bxes/ctn	S50021005
	20X20cm	5pcs/bx,40bxes/ctn	S50011006