

Alexer[®]

Polyurethane film backing

- Water proof, bacteria barrier
- Prevent external friction
- Managing the micro environment

Foam absorbent layer

- prevent wound maceration
- provide a moist wound environment
- reduce the risk of infection

Silicone contact layer


- Suitable for sensitive skin
- prevent adhere to the skin
- hydrophobicity which reduce the pain of dressing changes and avoid secondary damage

SILICONE FOAM DRESSING

◇ Without Border

DESCRIPTION

Alexer® silicone foam dressing-without border is useful for moderately to high exudating wounds, such as post-operative wounds, sacral bedsore wounds, pressure ulcers, diabetic foot ulcers, venous leg ulcers, neuropathic ulcers, arterial ulcers, post-operative, I - II degree burns etc.



Protective film


Breathable, waterproof, top film provides bacterial barrier.

Absorbent foam

Highly-absorbent, polyurethane foam core facilitates a wear time of up to 7 days.

Gentle adhesive

Easy to apply and remove, and shown to minimize pain and trauma during dressing changes.



Type	Size	Package	Ref.No.
 Square	10X10cm	5pcs/bx,40bxes/ctn	S10013001
	15X15cm	5pcs/bx,40bxes/ctn	S10013002
	5X5cm	5pcs/bx,40bxes/ctn	S10013003
	7.5X7.5cm	5pcs/bx,40bxes/ctn	S10013005
	10X20cm	5pcs/bx,40bxes/ctn	S10013006
	20X20cm	5pcs/bx,40bxes/ctn	S10013007