


Alexer[®]



SILICONE TAPE

APPLIED RANGE

Used to affix and fix dressings to wound surfaces or to fix other medical devices to specific parts of the human body.



Protect the heel



Fixed catheter



Scar repair



Correct sleep opening

PRODUCT FEATURES



- **MEDICAL-GRADE SILICONE** - Alexer® silicone tape is made with hospital medical-grade soft silicone. It's non-irritating and easily applied and removed, provides pain-free removal and reposition. It's also friendly to sensitive and fragile skin and maintains integrity during use, even when repositioned.
- **FLEXIBLE & EASY-TEAR** - No need for scissors to cut the tape- you can effortlessly apply the silicon.gel tape on to any part of your body. Furthermore, the flexible design makes for comfortable application to hard contours of the body.
- **CAN BE USED FOR MEDICAL APPLICATION** - Alexer® silicone tape is designed for the reduction of scars, which is ideal for a variety of scars, including surgical scars, c-section scars, and acne scars. Just place over old or new scars to reduce itching, redness, and discomfort. And also designed to be ideal for fixation of medical devices such as drains, tubes, probes, electrodes, and much more. It's designed to hold firmly and gently remove without pulling at the skin.

SPECIFICATIONS



Type	Size	Package	Ref.No.
Tape	2cmX1.5m	1 roll/box, 200rolls/ctn	S30011001
	4cmX1.5m	1 roll/box, 200rolls/ctn	S30011002
	4cmX3m	1 roll/box, 200rolls/ctn	S30011003
	2.5cmX3m	1 roll/box, 200rolls/ctn	S30011004

